

MacOS compatibility chart DinoXcope



	<i>Dino-Lite 311S/313/ 411/413/423 series</i>	<i>Dino-Lite 2xxx/3xxx series</i>	<i>Dino-Lite 3613/3713/ 3715 series</i>	<i>Dino-Lite 4013/4113/ 4023 series</i>	<i>Dino-Lite AM 4115/4515*/ 4815** series</i>	<i>Dino-Lite AF 4115*/4515*/ 4915**</i>	<i>AM4117/4517/ 4917</i>	<i>Dino-Lite 7013/7023 series</i>	<i>Dino-Lite 7115*/7515*/ 7815**/ 7915**/ 7025X**</i>	<i>Dino-Lite 73115*/73515*/ 73915** series</i>
Mac OS X 10.5 (Leopard)	Version 1.15 + Driver	Version 1.15	no	Version 1.15	no	no	no	Version 1.15	no	no
Mac OS X 10.6 (Snow Leopard)	Version 1.15 + Driver	Version 1.15	no	Version 1.15	no	no	no	Version 1.15	no	no
Mac OS X 10.7 (Lion)	Version 1.15 + Driver	Version 1.15	no	Version 1.15	no	no	no	Version 1.15	no	no
Mac OS X 10.8 (Mountain Lion)	Version 1.15 + Driver	Version 1.15	no	Version 1.15	no	no	no	Version 1.15	no	no
Mac OS X 10.9 (Mavericks)	Version 1.15 + Driver	Version 1.15	no	Version 1.15	Version 1.15	no	no	Version 1.15	no	no
Mac OS X 10.10 (Yosemite)	Version 1.15 + Driver	Version 1.15	no	Version 1.15	Version 1.15	no	no	Version 1.15	no	no
Mac OS X 10.11 (El Capitan)	no	Version 2.0.4	no	Version 2.0.4	Version 2.0.4	no	no	Version 2.0.4	Version 2.0.4	Version 2.0.4
Mac OS X 10.12 (Sierra)	no	Version 2.0.4	no	Version 2.0.4	Version 2.0.4	Version 2.0.4	no	Version 2.0.4	Version 2.0.4	Version 2.0.4
Mac OS X 10.13 (High Sierra)	no	Version 2.0.4	no	Version 2.0.4	Version 2.0.4	Version 2.0.4	no	Version 2.0.4	Version 2.0.4	Version 2.0.4
Mac OS X 10.14 (Mojave)	no	Version 2.1.3	no	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3
Mac OS X 10.15 (Catalina)	no	Version 2.1.3	no	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3
MacOS 11 (Big Sur)	no	Version 2.1.3	no	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3
MacOS 12 (Monterey)	no	Version 2.1.3	no	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3



	<i>Dino-Lite Medical DermaScopes MEDL4DW, MEDL7DW</i>	<i>Dino-Lite Medical DermaScopes MEDL4DM, MEDL7DM</i>	<i>Dino-Lite Medical EarScopes MEDL3E, MEDL4E, MEDL4EP</i>	<i>Dino-Lite Medical TrichoScopes, MEDL3H, MEDL4HM, MEDL4HMA, MEDL4HMA5, MEDL4HVW, MEDL7HM, MEDL7HMA, MEDL7HM4</i>	<i>Dino-Lite Medical CapillaryScopes, MEDL4N Pro*, MEDL4N5 Pro*</i>	<i>Dino-Lite Medical CapillaryScope MEDL4N5</i>	<i>Dino-Lite Medical IriScope MEDL4R</i>	<i>Dino-Lite Medical MEDL4FW</i>
Mac OS X 10.5 (Leopard)	no	Version 1.15	Version 1.15	Version 1.15	no	no	Version 1.15	Version 1.15
Mac OS X 10.6 (Snow Leopard)	no	Version 1.15	Version 1.15	Version 1.15	no	no	Version 1.15	Version 1.15
Mac OS X 10.7 (Lion)	no	Version 1.15	Version 1.15	Version 1.15	no	no	Version 1.15	Version 1.15
Mac OS X 10.8 (Mountain Lion)	no	Version 1.15	Version 1.15	Version 1.15	no	no	Version 1.15	Version 1.15
Mac OS X 10.9 (Mavericks)	Version 1.15	Version 1.15	Version 1.15	Version 1.15	Version 1.15	Version 1.15	Version 1.15	Version 1.15
Mac OS X 10.10 (Yosemite)	Version 1.15	Version 1.15	Version 1.15	Version 1.15	Version 1.15	Version 1.15	Version 1.15	Version 1.15
Mac OS X 10.11 (El Capitan)	Version 2.0.4	Version 2.0.4	Version 2.0.4	Version 2.0.4	Version 2.0.4	Version 2.0.4	Version 2.0.4	Version 2.0.4
Mac OS X 10.12 (Sierra)	Version 2.0.4	Version 2.0.4	Version 2.0.4	Version 2.0.4	Version 2.0.4	Version 2.0.4	Version 2.0.4	Version 2.0.4
Mac OS X 10.13 (High Sierra)	Version 2.0.4	Version 2.0.4	Version 2.0.4	Version 2.0.4	Version 2.0.4	Version 2.0.4	Version 2.0.4	Version 2.0.4
Mac OS X 10.14 (Mojave)	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3
Mac OS X 10.15 (Catalina)	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3
MacOS 11 (Big Sur)	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3
MacOS 12 (Monterey)	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3	Version 2.1.3

* FLC and AMR compatible from version 1.22 and up
 ** FLC and AMR compatible from version 1.22. No EDOF/EDR functions available